

SM-K-A ANESTHESIA TROLLEY

- Designed for convenient arrangement of instruments, medicines, and materials in medical facilities
- Made of premium stainless steel



SM-K-A
Dimensions (L/W/H):
- 1060x655x915 mm
Ø 125 mm castors



SM-K-A1
Dimensions (L/W/H):
- 760x655x915 mm
Ø 125 mm castors



LEARN MORE

SM-K MEDICAL TABLES

- Used for placing instruments, medical consumables, and medicines in dressing rooms, operating rooms, and other rooms of healthcare facilities
- Push bar for convenient transportation



SM-K-01
(stainless steel)
SM-K-11
(steel with polymeric coating)
Dimensions (L/W/H):
- 720x460x805 mm
Ø 80 mm castors



SM-K-02
(stainless steel)
SM-K-12
(steel with polymeric coating)
Dimensions (L/W/H):
- 720x460x805 mm
Ø 80 mm castors



SI MEDICAL INSTRUMENT TABLES

- Designed to equip surgical departments of healthcare facilities
- Made of premium stainless steel



SI-51
Dimensions (L/W/H):
- 700x500x900 mm
Adjustment by foot pedal
Height adjustment range:
900...1300 mm
Ø 80 mm castors



SI-50
Dimensions (L/W/H):
- 700x500x950 mm
Manual adjustment
Height adjustment range:
900...1250 mm
Ø 80 mm castors



SI-60
Dimensions (L/W/H):
1000x650x850 mm



SI-60-1
Dimensions (L/W/H):
1000x650x850 mm
(without castors)



SI-61
Dimensions (L/W/H):
1700x650x850 mm

SI-60-2
Dimensions (L/W/H):
1200x675x850 mm
Ø 80 mm castors



SI-62
Dimensions (L/W/H):
1200x725x1000 mm
Ø 150 mm castors

SI-61-1
Dimensions (L/W/H):
1700x650x850 mm
(without castors)

SI-61-2
Dimensions (L/W/H):
1800x675x850 mm
Ø 80 mm castors



SM-K-03
(stainless steel)
SM-K-13
(steel with polymeric coating)
Dimensions (L/W/H):
- 720x460x905 mm
Ø 80 mm castors



SM-K-04
(stainless steel)
SM-K-14
(steel with polymeric coating)
Dimensions (L/W/H):
- 720x460x905 mm
Ø 80 mm castors



SM-K-05
(stainless steel)
SM-K-15
(steel with polymeric coating)
Dimensions (L/W/H):
- 1040x700x805 mm
Ø 100 mm castors



SM-K-06
(stainless steel)
SM-K-16
(steel with polymeric coating)
Dimensions (L/W/H):
- 1040x700x805 mm
Ø 100 mm castors



SM-K-07
(stainless steel)
SM-K-17
(steel with polymeric coating)
Dimensions (L/W/H):
- 1040x700x905 mm
Ø 100 mm castors



SM-K-08
(stainless steel)
SM-K-18
(steel with polymeric coating)
Dimensions (L/W/H):
- 1040x700x905 mm
Ø 100 mm castors