



SMP MEDICAL STOOLS



- Stool base is made of premium stainless steel/aluminum
- Molded seats and backrests are made of integral polyurethane foam
- Color of seat and backrest covers may be chosen (accessory)



SMPS

Ø370 mm seat
Seat height adjustment range:
430...630 mm
Screw-type mechanism
Weight max. 12 kg
Ø 50 mm castors



SMP

Ø370 mm seat
Seat height adjustment range:
430...630 mm
Screw-type mechanism
Weight max. 7,5 kg
Ø 50 mm castors



SMPP

Ø370 mm seat
Seat height adjustment range:
490...630 mm
Hydraulic jack
Weight max. 14 kg
Ø 50 mm castors



SMPP-01

Ø370 mm seat
Seat height adjustment range:
490...630 mm
Hydraulic jack
Weight max. 12 kg
Ø 50 mm castors



SMPP-02

Ø370 mm seat
Seat height adjustment range:
480...670 mm
Pneumatic drive
Weight max. 13 kg
Ø 50 mm castors



SMPP-03

Ø370 mm seat
Seat height adjustment range:
480...670 mm
Pneumatic drive
Weight max. 11 kg
Ø 50 mm castors

Cover colors

