



NM

MEDICAL NEGATOSCOPES



LEARN MORE

- Antistatic anti-glare glass
- Table or wall mounted
- Independent smooth brightness adjustment for each screening area
- Independent switching on of each screening area
- Low power consumption



Dimensions	HM-1	HM-2	HM-3
Length, mm	405	765	1125
Width, mm	115	115	115
Height, mm	515	515	515
Weight, kg	8	13	18

Max. brightness: (4,200...1,000) cd/m²
Single screening area size: 360x400 mm
Start-up time: 5 s