





- Used to house medical diagnostic equipment, instruments, tools, and other medical devices
- Made of engineering steel and aluminum profile coated with polymer paint
- Strong, durable construction, convenient for sanitary treatment

SA-1
Dimensions (L/W/H):
- 700x660x930 mm
Shelf dimensions (W/D):
- 550x490 mm



Accessories



SA-2
Dimensions (L/W/H):
- 700x660x1200 mm
Shelf dimensions (W/D):
- 550x490 mm



SA-3
Dimensions (L/W/H):
- 700x660x1470 mm
Upper shelf dimensions (W/D):
- 550x300 mm
Shelf dimensions (W/D):
550x490 mm



SA-4
Dimensions (L/W/H):
- 550x655x910 mm
Shelf dimensions (W/D):
- 490x400 mm



SA-5
Dimensions (L/W/H):
- 550x655x1180 mm
Shelf dimensions (W/D):
- 490x400 mm



SA-8
Dimensions (L/W/H):
- 1000x660x1200 mm
Shelf dimensions (W/D):
- 850x490 mm



SA-9
Dimensions (L/W/H):
- 1000x660x1470 mm
Upper shelf dimensions (W/D):
- 850x300 mm
Shelf dimensions (W/D):
- 850x490 mm

