


MEDIN-ELPIS


RESUSCITATION STRETCHER



- Comfortable and maneuverable patient transporter with foot-operated hydraulic height and longitudinal tilt adjustment
- May be operated from both sides
- X-ray transparent table top with X-ray cassette tunnel
- Removable mattress dimensions: 1,950x640x80 mm



Accessories




MEDIN-ELPIS 02

RESUSCITATION STRETCHER



- Comfortable and maneuverable patient transporter with foot-operated hydraulic height adjustment
- May be operated from both sides
- X-ray transparent table top with X-ray cassette tunnel
- Removable mattress dimensions: 1,950x640x80 mm



Accessories

